

My Quality Solutions




To generate confidence



Goal



Materials



Reference



Frequency



PRECISION

- Evaluation of the analytical performance and consistency of a testing system (Reagents + Instruments)

- 2 or 3 levels of assayed QC
- Similar to patient samples

- Self evaluation
- Manufacturer ranges or lab-established ranges

- Daily



RELATIVE PRECISION

- Relative performance and bias of a lab's results compared to peers

- 2 or 3 levels of assayed IQC

- Peer group
- Self evaluation

- Real time



ACCURACY

- Independent and blind verification of analytical performance of a laboratory

- Multiple non assayed levels
- Blind samples
- Similar to patient samples

- Peer group

- Periodically

Features	My Daily QC	My Expert QC	My Qualiris QC
Contribution to Lab QC Management	Stago IQC Management	Externalised IQC	External Quality Assessment
Web based solution		✓	✓
Automated upload of QC results		✓	
Easy interpretation (alerts and flags)	✓	✓	✓
QC range management	✓	✓	
Meaningful graphs and statistical indicators	✓	✓	✓
QC Lot management	✓	✓	
Automated management of data	✓	✓*	
Automated generation of report		✓	✓
Contribution to lab accreditation (ISO 15189)	✓	✓	✓
Expert support	✓	✓	✓
Comprehensive offer (routine and special assays)	✓	✓	✓
Real time peer group statistics		✓	
Peer group data comparison		✓	✓
Multiple lab management		✓	✓

* On Max Generation analysers. Other data input methods available for other Stago analysers.

Connect.One, Qualiris, My Expert QC, STA, STA Coag Expert are Stago Group trademarks. The rights concerning the tradenames and logos used in this document are held by the Stago Group. These trademarks may not be used without authorisation of the Stago Group.

Conception: L2R.fr ©2017 Diagnostica Stago - All rights reserved - Non-contractual photos - 11/2017 - Ref. 300361
This document contains information on products which is targeted to a wide range of audiences and could contain product details or information otherwise not accessible or valid in your country.

Increase confidence Compare results with global peers

For further information, please contact:



At the Heart of Haemostasis

Diagnostica Stago S.A.S.
RCS Nanterre B305 151 409
3, allée Thérèse
92600 Asnières sur Seine
France
Ph.: +33 (0)1 46 88 20 20
Fax: +33 (0)1 47 91 08 91
webmaster@stago.com
www.stago.com



My Expert QC: the Externalised IQC offer by Stago

My Expert QC is a real time external quality control program that continuously monitors the performance of your analyser based on the routine use of Stago internal quality control materials.

My Expert QC offers a flexible group comparison from local to global, covering all Stago QC from routine to specialty.

How to export my QC data?

- My Expert QC operates via Connect.One, the new communication platform dedicated to Haemostasis
- Daily QC data are transferred by secure internet connection directly to My Expert QC Servers

How can My Expert QC help me?

Confident results

- Supports your lab to meet accreditation & compliance requirements
- Daily IQC monitoring provides higher level of confidence in system performance
- Extended comparison groups (Worldwide and Nationwide) for higher statistical power of system trueness
- Additional levels of comparison (global, country or laboratory network) for a deeper analysis

Stress-free QC management

- Expert support from the Stago team to assist in interpretation and troubleshooting
- Robust peer group comparison statistics with true value-added indicators
- Meaningful reports and graphs to complete your quality management approach

Cost effective solution

- Enable proactive intervention as soon as trends are observed
- Group data are generated real-time to reduce unnecessary troubleshooting


Easy QC management

- Simple and intuitive interface
- Simple reports with relevant statistical analysis
- On-demand, customised group comparison reports
- Access from any computer connected to internet



How to access my QC data?

- Access from any computer connected to the internet
- Password-protected, laboratory-specific data are available on: www.myexpertqc.com
- Quality Managers can view results from their networked analysers with a single login



What types of report are available?

- **Quick report**
At-a-glance report compares your laboratory with the peer group
- **Full report**
Full report provides additional statistics and comparison levels (ex. comparison with all STA analysers, worldwide comparison and hospital systems comparison)
- **Graphical report**
Displays the graphical distribution of results for the group and the relative position of your laboratory or analyser within the group